



Recombinant Mouse Peroxiredoxin-4 (Prdx4)

Product Code	CSB-EP522510MO-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O08807
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	LESDERFRTR ENECHFYAGG QVYPGEASRV SVADHSLHLS KAKISKPAPY WEGTAVINGE FKELKLTIDYR GKYLVEFFYP LDFTFVCPTE IIAFGDRIEE FKSINTEVVA CSVDSQFTHL AWINTPRRQG GLGPIRIPLL SDLNHQISKD YGVYLED SGH TLRGLFIIDD KGVLRQITLN DLPVGRSVDE TLRLVQAFQY TDKHGEVCPA GWKPGSETII PDPAGKLYF DKLN
Source	E.coli
Target Names	Prdx4
Protein Names	Recommended name: Peroxiredoxin-4 EC= 1.11.1.15 Alternative name(s): Antioxidant enzyme AOE372 Peroxiredoxin IV Short name= Prx-IV Thioredoxin peroxidase AO372 Thioredoxin-dependent peroxide reductase A0372
Expression Region	41-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is an antioxidant enzyme and belongs to the peroxiredoxin family. The protein is localized to the cytoplasm. Peroxidases of the peroxiredoxin family reduce hydrogen peroxide and alkyl hydroperoxides to water and alcohol with the use of reducing equivalents derived from thiol-containing donor molecules. This protein has been found to play a regulatory role in the activation of the transcription factor NF-kappaB.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.