



Recombinant Rat Intestinal-type alkaline phosphatase 1 (Alpi)

Product Code	CSB-BP001627RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15693
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	VIPVEEENPV FWNQKAKEAL DVAKKLQPIQ TSAKNLILFL GDGMGVPTVT ATRILKGQLG GHLGPETPLA MDHFPFTALS KTYNVDRQVP DSAGTATAYL CGVKANYKTI GVSAAARFNQ CNSTFGNEVF SVMHRAKKAG KSVGVTTR VQHASPAGTY AHTVNRDWYS DADMPSSALQ EGCKDIATQL ISNMDIDVIL GGGRKFMFPK GTPDPEYPGD SDQSGVRLDS RNLVEEWLAK YQGTRYVWNR EQLMQASQDP AVTRLMGLFE PTEMKYDVNR NASADPSLAE MTEVAVRLLS RNPQGFYLFV EGGRIDQGHH AGTAYLALTE AVMFDSAIEK ASQLTNEKDT LTLITADHSH VFAFGGYTLR GTSIFGLAPL NAQDGKSYTS ILYGNGPGYV LNSGNRPNVT DAESGDVNYK QQAAVPLSSE THGGEDVAIF ARGPQAHLVH GVQEQNYIAH VMAFAGCLEP YTDCGLAPPA DENRPTTPVQ N
Source	Baculovirus
Target Names	Alpi
Protein Names	Recommended name: Intestinal-type alkaline phosphatase 1 Short name= IAP-1 Short name= Intestinal alkaline phosphatase 1 EC= 3.1.3.1 Alternative name(s): Intestinal alkaline phosphatase I Short name= IAP-I
Expression Region	21-511
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	There are at least four distinct but related alkaline phosphatases: intestinal, placental, placental-like, and liver/bone/kidney (tissue non-specific). The intestinal alkaline phosphatase gene encodes a digestive brush-border enzyme. This enzyme is upregulated during small intestinal epithelial cell differentiation.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.