



# Recombinant Schizosaccharomyces pombe 1,3-beta-glucanosyltransferase gas5 (gas5)

<b>Product Code</b>	CSB-MP522576SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O13692
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	DSASSAIK IKGNAFFNSD TNERFYVRGV DYQPGGSSTL VDPLADTSIC KRDLPYLQGL NINTIRVYQV DNSANHDECM SALQDAGIYV ILDLATSSNS ISRLDAASSY NAVFLQGIFA TIDAFKNYTN VLGFFAGNEV ANTAENSATT TWVKAALRDA KEYISKNSDR DIPVGYSAAD VAEIRVQCAD FFACGNSSVR ADFYGMNMYE WCGADSSFTI SGYDQRMEEF ANYSIPLFLS EYGCNDVTKE SDGTPDRPFD EVDAIFSSEM SSVFSGGLVY QYSEEGNNYG LVVIDGDNVT ISKNYETLKE KYASAANYTG DGDYSSSPAT LTCPADDSYF TSFPLPTMPS EAKGFIESGA GQPLGFNAPS NQEFSANATA LVSPGPHSVS TTINTNIVQA TISQSSTSGS SSGSSSASTT ASSSSVSSGS SISSGSSSMS TSYTSASGSS AHSSGSSSGS SSATS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	gas5
<b>Protein Names</b>	Recommended name: 1,3-beta-glucanosyltransferase gas5
<b>Expression Region</b>	23-485
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.