



Recombinant Schizosaccharomyces pombe ATP-dependent RNA helicase uap56 (uap56)

Product Code	CSB-EP522590SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O13792
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MASAQEDLID YEEEEELVQD QPAQEITPAA DTAENGEKSD KKGSYVGIHS TGFRDFLLKP ELLRAITDSG FEHPSEVQQV CIPQSILGTD VLCQAKSGMG KTAVFVLSTL QQIEPVDGEV SVLVLCHTRE LAFQIKNEYA RFSKYLDPVR TAVFYGGINI KQDMEAFKDK SKSPHIVVAT PGRLNALVRE KILKVN SVKH FVLDECDKLL ESVDMRRDIQ EVFRATPPQK QVMMFSATLS NEIRPICKKF MQNPLEIYVD DETKLT LHGL QQHVKLEEK AKNRKINDLL DSLEFNQVVI FVKSVSRANE LDROLLRECNF PSICIHGGLP QEERIKRYKA FKDFDKRICV ATDVFGRGID IERVNIVINY DMPDSPDSYL HRVGRAGRFG TKGLAITFSS SEEDSQILDK IQERFEVNIT ELPDEIDVGS YMNA
Source	E.coli
Target Names	uap56
Protein Names	Recommended name: ATP-dependent RNA helicase uap56 EC= 3.6.4.13
Expression Region	1-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.