



Recombinant Schizosaccharomyces pombe Serine/threonine-protein kinase sid1 (sid1)

Product Code	CSB-EP522670SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O14305
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MYPLNANSYT LLRKLGSGSF GVVWKARENV SGDIIAIKQI DLETGIDDIT DIEQEVFMLS NCNSSNVIQY YGCFVDGYTL WILMEHMDGG SVSGLLKMGR LNEQVISIIL REVLYGLNYL HGQNKIHRDI KAA NILSSS TGNVKLADFG VAAQLSNAAS RRHTFVGTPF WMAPEVIQQT SYGLAADIWS LGITAIEMAN GIPPRATMHP MRVIFEIPQS EPPKLDDHFS PTFRDFVSCC LDLNPNMRWS AKELLQHPFI KSAGTVKDII PLLVQKENKL FDDSDQSVLE ETINNTLKPF EEPIAEGNAD IEDWTFETVK KSDSTVLGNT SIPKNSIISS QNKEELPSSI KYLEKTIMSD QATPHPFSSK LSEKSSSYHK SLTSDFAMKH YIKSTIRSML LNDKLSATQR SSLESFYTSF ISL DKNLSSK FVNQITPDNR LHHKKQKRSP ISQLLFSRWL EETEKRRSLN G
Source	E.coli
Target Names	sid1
Protein Names	Recommended name: Serine/threonine-protein kinase sid1 EC= 2.7.11.1 Alternative name(s): STE20-like kinase sid1
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.