



Recombinant Uncharacterized protein F52F12.7 (F52F12.7)

Product Code	CSB-YP522718CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O17883
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	EKTNPNLTLH DVTDKLEGDD EKYATALKTC GEVFAEVEAI FNDENYLSHA GWFKDESNNE GDVVYAKDTP HGRMVTISTE LPMPVEDVMK ETWNGMEALP EWNQNINFAA RIAAPTSNFD IVTYGNNDVL VVSGREFVSA RIWRKVG DGF ILASRSVTVP SFSKSHKGKV RAHLHLAGAR FRPNPENPET TLTDVVMLAD LKGYLPKMIV NQVIGRIMIM DTVTNRRHFQ NLKAKRTNN
Source	Yeast
Target Names	strl-1
Protein Names	Recommended name: Uncharacterized protein F52F12.7
Expression Region	21-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.