



Recombinant Nuclear hormone receptor family member nhr-79 (nhr-79)

Product Code	CSB-YP522722CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O18141
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MVRGKCMVCD SPNATNYHFG AQSCKACAAF FRRSIAMDQC YECLGDGVHS CKIDHTLRN CRHCRLKKCQ KAGMMRDLVQ AKREIKSDKG KNSRNSSQSE DFFSPPPEQP GPSNYFDQFP WQASENEMDP YPSSPKIRKM SENLMGIDEL QQFIQIPAGV DDFPDIFTDS DEIRSRMTSI SEASAYVTGM EEEEERLFLA GLYTEQVINL NMRRRITYTD RLLGSVFDAP CVCPYDKADL KLFDRHTYRQ KNRNDYTMIL DYINRFPEFE SLSKSEKTVL FRATAAVDVL LDQSYYSQVI FPTEDVLVTA NGEYLPMPNPM PKVENQRDSG NFHSEDEDYDR FKMLTSMKVR QWLHVCEPMK KLDMSLAEFS LFKALTIWHY NYYKLQCTGK QICSRQRDDI FRTLLICDD EGHDSALIRA SEIVLAVGIA MAEVHEMVT YIEITVYDVL DDPILKDMMLK FQY
Source	Yeast
Target Names	nhr-79
Protein Names	Recommended name: Nuclear hormone receptor family member nhr-79
Expression Region	1-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.