



Recombinant Arabidopsis thaliana Putative serine carboxypeptidase-like 30 (SCPL30)

| | |
|--------------------------|--|
| Product Code | CSB-BP522808DOA |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | O23364 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | CA RQHRRSFLAK EADLVTNLPG QPDVSFKHYA GYVPVDKSNG RALFYWFFEA MDLPKEKPLV LWLNNGPGCS SVGYGATQEI GPFLADTNEK GLIFNPYAWN KEVNMLFLES PVGVGFSYSN TSSDYLNLDH HFAKKDAYTF LCNWFEEKFPE HKGNEFYIAG ESYAGIYVPE LAELVYDNNE KNNDLSLHIN LKGFLGPNPD ISNPDDWRGW VDYAWSHAVI SDETHRNINR LCNFSSDDVW NNDKCNEAIA EVDKQYNEID IYSLYTSACK GDSAKSSYFA SAQFKTNYHI SSKRMPPRRL AGYDPCLEDDY VKVYYNRADV QKALHASDGV NLKNWSICNM EIFHNWTYVV QSVLPIYQKL IAGGLRIWVY SGTDTGCIPV LGTRYSLNAL GLPIKTAWRP WYHEKQVSGW VQEYDGLTFA TFRGAGHTVP SFKPSSSLAF ISAFVKGVPPL SSSRVETN |
| Source | Baculovirus |
| Target Names | SCPL30 |
| Protein Names | Recommended name: Putative serine carboxypeptidase-like 30 EC= 3.4.16.- |
| Expression Region | 29-488 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |