



Recombinant Arabidopsis thaliana Putative serine carboxypeptidase-like 30 (SCPL30)

Product Code	CSB-MP522808DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O23364
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	CA RQHRRSFLAK EADLVTNLPG QPDVSFKHYA GYVPVDKSNG RALFYWFFEA MDLPKEKPLV LWLNGGPGCS SVGYGATQEI GPFLADTNEK GLIFNPYAWN KEVNMLFLES PVGVGFSYSN TSSDYLNLDL HFAKKDAYTF LCNWFEEKFPE HKGNEFYIAG ESYAGIYVPE LAELVYDNNE KNNDLSLHIN LKGFLGPNPD ISNPDDWRGW VDYAWSHAVI SDETHRNINR LCNFSSDDVW NNDKCNEAIA EVDKQYNEID IYSLYTSACK GDSAKSSYFA SAQFKTNYHI SSKRMPPRRL AGYDPCLEDDY VKVYYNRADV QKALHASDGV NLKNWSICNM EIFHNWTYVV QSVLPIYQKL IAGGLRIWVY SGDTDGCIPV LGTRYSLNAL GLPIKTAWRP WYHEKQVSGW VQEYDGLTFA TFRGAGHTVP SFKPSSSLAF ISAFVKGVPPL SSSRVETN
Source	Mammalian cell
Target Names	SCPL30
Protein Names	Recommended name: Putative serine carboxypeptidase-like 30 EC= 3.4.16.-
Expression Region	29-488
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.