



Recombinant Arabidopsis thaliana Peroxidase 41 (PER41)

Product Code	CSB-BP522815DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O23609
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	APPPNLT KDY YQKTCPDFNK IVRETVTPKQ GQQPTTAAGT LRLFFHDCFM EGCDASVLIA TNSFNKAERD DDLNESLPGD AFDIVTRIKT ALELSCPGVV SCADILAQAT RDLVTMVGGP FYEVKLGRKD GFESKAHKVK GNLPLANQSV PDMLSIFKKN GFTLKELVAL SGGHTIGFSH CKEFSNRIFP KVDPELNAKF AGVLKDLCKN FETNKTMAAF LDPVTPGKFD NMYFKNLKRG LGLLASDHIL FKDPSTRPFV ELYANNQTAF FEDFARAMEK LGRVGVKGEK DGEVRRRCDH FNKLVN
Source	Baculovirus
Target Names	PER41
Protein Names	Recommended name: Peroxidase 41 Short name= Atperox P41 EC= 1.11.1.7
Expression Region	21-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.