



Recombinant Arabidopsis thaliana Ras-related protein RABE1a (RABE1A)

Product Code	CSB-YP522840DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O24466
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAAPPARARA DYDYLIKLLL IGDSGVGKSC LLLRFSDGSF TTSTFITTIGI DFKIRTIELD GKRIKLQIWD TAGQERFRTI TTAYYRGAMG ILLVYDVTDE SSFNNIRNWI RNIEQHASDS VNKILVGNKA DMDESKRAVP KSKGQALADE YGMKFFETSA KTNLNVVEEVF FSIKDIKQR LADTDARAEP QTIKINQSDQ GAGTSQATQK SACCCT
Source	Yeast
Target Names	RABE1A
Protein Names	Recommended name: Ras-related protein RABE1a Short name= AtRABE1a Alternative name(s): Ras-related protein Rab8B Short name= AtRab8B
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.