



Recombinant Schizosaccharomyces pombe Uncharacterized but2-like protein C27D7.09c (SPAC27D7.09c, SPAC27D7.10c)

Product Code	CSB-BP523649SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O42663
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	SPLHKRD DPLNNKQFGL MSIHSGNSYV HLHPFYVGDS GAIYLDPTDG TSTSAVFSMS NGRLVVGPLY ASVDSNGTTI FKSDANAASST KFSAGDATNV GYNLLYNGTQ SAVACPASDN DQVYQVYFGA GNGNPNCAGI AIEAFVSPSS SSSSSSSAAT STSTRVSSSA KASTSSGAIA YTTKCVVVPV TASATATAKA ASAAASSAVY PLFPHGIRLI DSANPSSNSG NVYSPVVFQK QNNHTNTIFT FDVPQVSGSC ELNLHLDTSR FPITVEGSNG VGQFILFNLS SVANDSTVYS NRPNRIAEIG RFNCSSSGCD YATNVTC PNS YTAVSYEMMA LTDDSYLSFF EEADPLEGLT LRV
Source	Baculovirus
Target Names	SPAC27D7.09c
Protein Names	Recommended name: Uncharacterized but2-like protein C27D7.09c
Expression Region	24-383
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.