



Recombinant Schizosaccharomyces pombe RNA-binding post-transcriptional regulator cip1 (cip1)

| | |
|--------------------------|---|
| Product Code | CSB-EP523666SXV-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | O42923 |
| Product Type | Recombinant Protein |
| Immunogen Species | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MPNEIFPWKI RVDESRGLAG NSGTKKSNHE ALSRLQSPLN SPKLQPIGSP QASRKTSGSG SSAPLYPKWS GALSLASSRA ASPAPSDSFP TFGYSQLGGL ENSSKGSALF NSANSIGTPY LSSRNSNSAN EASAMAFHNV SPPSGAESSS ESKSFSASGK GNKADTSAEP SLDAFNSTQI KAGSTANSNS TPVEPGEDTI PTAIVVKNIP FSLEKDTLLD HFKQLGIPRP YAFNYHYDNG IFRGLAFANF YRPEEAQVVV QTLNGYEING RRLRVEWKRQ LPAAERЕКVE KAKKRQAEER RRKQQYKMFЕ VSFTDQGLNL NDTGTLETYS RLLLFАHQCI PSREITFETT SKDGNLLNAI RIFCLYFDLD YYARPNGEVL KLVVTHPNKK NTSVSQSQPA SPNLRFNMPA PLATRFLQEH SLNGTKSAPI TPPPSFAVPL TNQLRSIDDK IYGNESPLQK ASTLSSPFNS KNDNDASTSA SKQSFGVSHF |
| Source | E.coli |
| Target Names | cip1 |
| Protein Names | Recommended name: RNA-binding post-transcriptional regulator cip1 Alternative name(s): Csx1-interacting protein 1 |
| Expression Region | 1-490 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |