



Recombinant *Schizosaccharomyces pombe* Splicing factor U2AF-associated protein 2 (uap2)

Product Code	CSB-EP523681SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O43120
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSSQPFWDER IHRWRCLGSE GNELIYIDEE QTWKDYDPNS LKMNKAGSTG AEVSDVTAEA TEGKESSNGE DRHTKRLYES TSAEGYPSGS RNKKSSENS EASPAPVINK AVYIQGLPLD VTVDEIEEVF KKCGVIAKNI DNGTPRIKIY RTEDGTPKGD ALIVFFRSES VELAEQLFDD TEFRYGSGQK MRVQKANIDY KKEKTVNKDV GGALKKKALR LRQQMQQIS SWDDVDEEVD DKRKKRFNKI VVLKHIFTLE ELDKTPELLI DLKDDITEEA EKCGRVTVNVV LYDKEPDGVV TVRFSNNEEA EACVRLMQGR YFDGRVVEAS IYDGKVRFQK SGKHTLDDEE DEEKRLEKFA DWLENSN
Source	E.coli
Target Names	uap2
Protein Names	Recommended name: Splicing factor U2AF-associated protein 2 Alternative name(s): Cold sensitive U2 snRNA suppressor 2 homolog
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.