



Recombinant Arabidopsis thaliana Putative F-box/kelch-repeat protein At4g35120 (At4g35120)

Product Code	CSB-EP523736DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O49618
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MDTAKVEPPQ ETKKTILTAQ TPLSMSISL PDEIVLSFLA LISKSYYRSL SLVSKSFYSL LSSTEIYAAR SHIGATEPCP YVCLWLPPKH RWFTLAEIEG KLSLEPVRLS SSYPRTVNS TTVAAGTEIY KIGGTVKGKR SRAVFLDCW THRWRAPNM RVS RVGAKSC FLDGNIYVIG GCRKSEESM NCGEVFDLKT QTWNPLSPS VNYAVHSNHK VAVSGERLYV ITRNNYAYD PNEGRWLDPV GSVDLQPITG PWSGGIEKVM KPITGRPWKY TWYSSSHRAW QRVMGLDVLV DKRGCGYRTI QLVNYGGKLL IIVSEWIMIL DGYSLIRSQE KEIWCAVIRL EERMSYFGPQ IWGEVESC NVVPSVPKSYQ LSSCQCVSV
Source	E.coli
Target Names	At4g35120
Protein Names	Recommended name: Putative F-box/kelch-repeat protein At4g35120
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.