



# Recombinant Pyrococcus horikoshii Peptide chain release factor subunit 1 (prf1)

<b>Product Code</b>	CSB-MP523929FHX
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O59264
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pyrococcus horikoshii (strain ATCC 700860 / DSM 12428 / JCM 9974 / NBRC 100139 / OT-3)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSKMTRHDAQ LYELKKKIEE LKKIRGRGTE LISLYIPAGY DLSKVMQQLR EEYSTAQNIK SKTTRKNVLG ALERAMQHLK LYKQTPENGL ALFVGNVSEI EGNTDIRLWA IIPPEPLNVR LYRCDQTFVT EPLEMLRVK DAYGLITVEK NEATIGLLRG KRIELDELDT SNVPGKTRAG GQSARRYERI REQETHEFMK RIGEHANRAF LPLLEKGELK GIIIGGPGPT KEEFVEGEYL HHELKKKILG VVDISYHGEY GLRELVEKAS DILRDHEVIR EKKLVNEFLK HVVKDTGLAT YGEREVRRAL EIGAVDVLLI SEGYNKVRVR VKCNHCGWEE LKTMSEEEYE VYRKKITKCP KCGSQNLTIK KWDVAEELIK MAEEAGSDVE IISLDTEEGQ QFYRAFGGLG AILRFKI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	prf1
<b>Protein Names</b>	Recommended name: Peptide chain release factor subunit 1 Alternative name(s): Translation termination factor aRF1
<b>Expression Region</b>	1-417
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.