



Recombinant *Pyrococcus horikoshii* Signal recognition 54 kDa protein (srp54)

Product Code	CSB-YP523932FHX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O59307
Product Type	Recombinant Protein
Immunogen Species	<i>Pyrococcus horikoshii</i> (strain ATCC 700860 / DSM 12428 / JCM 9974 / NBRC 100139 / OT-3)
Purity	≥85% (SDS-PAGE)
Sequence	MVLDSLGRAL SNALKKIARA GSVDEALVKE VVRDIQRALL QADVNVRLVL KLTKEIQRRA LEEKPPAGIS KKEHIIKIVY EELTKFLGTE AKPIEIKEKP TVLLTVGVQG SGKTTTVAKL ARYFQKRGYK VGVVCSDTWR PGAYHQLKQL LDPYHIEVFG DPNEKDAVKL AKEGVEYFKT RDVDLIIVDT AGRHKEEKDL IEEMRMISEE IRPHEVILVI DGTIGQQAYN QALAFKEATP IGSIIIVTKLD SSAKGGGALS AVAATGAPIK FIGVGEKIDD LEPFDPARFV SROLLGLGDIQ GLLEKFKELE KEVEFTEEDI DRFLRGKFTL KDMYAQLEAM RKMGPLKQIL RMIPGLGYSL PDELISVGEE RLRKFKVIMD SMTEEELMNP EIINYSRIKR IARGSGTSIK DVKELLTQYN QMKKFFKSMN KRQLSRLARR FGM
Source	Yeast
Target Names	srp54
Protein Names	Recommended name: Signal recognition 54 kDa protein Alternative name(s): SRP54
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.