



Recombinant Pyrococcus horikoshii 50S ribosomal protein L22P (rpl22p)

Product Code	CSB-MP523937FHX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O59423
Product Type	Recombinant Protein
Immunogen Species	Pyrococcus horikoshii (strain ATCC 700860 / DSM 12428 / JCM 9974 / NBRC 100139 / OT-3)
Purity	>85% (SDS-PAGE)
Sequence	MGKRFGYSFQ NFDPERMARA SGRDLRISPK LAVEVCRELR GMMLNDALRY LDEVIALKRP VPLRRYND SQ GHKPGKGF GP GRYPVKVAKA IKKVLLNAKN NAEQKGLDPD KLKIIHIAAH RGPVLRGWYP RAFGRATPFN EQTTHIEVVV EEIRR
Source	Mammalian cell
Target Names	rpl22
Protein Names	Recommended name: 50S ribosomal protein L22P
Expression Region	1-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.