



# Recombinant Schizosaccharomyces pombe U3 small nucleolar RNA-associated protein 6 (utp6)

<b>Product Code</b>	CSB-YP523984SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O60188
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAEKVQYYME QSVPELEDLL EKNIFNRDEI NNIKTRRVF EEKLARRQVK LNDFLSYIQY EINLETLRAK RHKRLNITGK ITISDYAGPR KVLFLFLRAT NKFFGDVTLW LDYIHYAQKI KAVNIVGKIC VAALQKHPNN AELWVVACDH EFSINANVSA ARALMNRALR LNQENPVIWA AYFRLELSYM TKLFARSQIL TGNISSKTEN ITNGVSEDTI GSLSSDTIQL PMVSMEEFLG SSSSEVRKND SDLNISDDIG NISSKEQQTQ KFANVLLQII LNSRKNLSLQ NYVGGFFVSVL DALFECFDVP VVQYMYQENI IGICNEHFEH FKNESGEIYG VLLHRWCFLE IFIKLRGAYP SNDSGISGKG IFGLKKN DPR ITSMVLTPG FVDDLQAVE KYQSISSDFQ IPFKTKKIFY SFFVKTLHAI SSSSAAESSI ALALHMLVIN AFQSMEKTKI LTFDDSLQEI YKEAQIQSGT FMSNATLS
<b>Source</b>	Yeast
<b>Target Names</b>	utp6
<b>Protein Names</b>	Recommended name: U3 small nucleolar RNA-associated protein 6 Short name= U3 snoRNA-associated protein 6
<b>Expression Region</b>	1-488
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.