



# Recombinant Human ProSAP-interacting protein 1 (PROSAPIP1)

<b>Product Code</b>	CSB-EP523986HU-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O60299
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAKLETLPVR ADPGRDPLLA FAPRPSELGP PDPRLAMGSV GSGVAHAQEF AMKSVGTRTG GGGSQGSFPG PRGSGSGASR ERPGRYPSED KGLANSLYLN GELRGSDHTD VCGNVVGSSG GSSSSGGSDK APPQYREPSH PPKLLATSGK LDQCSEPLVR PSAFKPVVVK NFHSMQNLCP PQTNGTPEGR QGPGGLKGGL DKSRTMTPAG GSGSGLSDSG RNSLTSLPTY SSSYSQHLAP LSASTSHINR IGTASYGSGS GGSSGGGSGY QDLGTSDSGR ASSKSGSSSS MGRPGHLGSG EGGGGGLPFA ACSPSPSAL IQELEERLWE KEQEVAALRR SLEQSEAAVA QVLEERQKAW ERELAE LRQG CSGKLQQVAR RAQRAQQGLQ LQVLR LQQDK KQLQEEAARL MRQREELEDK VAACQKEQAD FLPRIETKW EVCQKAGEIS LLKQQLKDSQ ADVSQKLSEI VGLRSQLREG RASLREKEEQ LLSLRDSFSS KQASLELGEG ELPAACKPA LTPVDPAEPQ DALATCESDE AKMRRQAGVA AAASLVSDG EAEAGGESGT RALRREVGRL QAELAAERRA RERQGASFAE ERRVWLEEKE KVIEYQKQLQ LSYVEMYQRN QQLERLRER GAAGGASTPT PQHGEKKAW TPSRLERIES TEI
<b>Source</b>	E.coli
<b>Target Names</b>	LZTS3
<b>Protein Names</b>	Recommended name: ProSAP-interacting protein 1 Short name= ProSAPiP1
<b>Expression Region</b>	1-673
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .