



Recombinant Arabidopsis thaliana FBD-associated F-box protein At1g61320 (At1g61320)

Product Code	CSB-BP524042DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O64788
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MAPPTKRTRV EMAESSNKRM KPSETVPEDV LELMMSTYLP VQSLLTTRVL SKRFRETEVR SLDLDFSGIY SRRRRKVEVV GIIEKVFNQH KGSEIHRFVL LLNHIGVEDK IISWTNTCFD KNIKELVLDV SKSRKVMAIP IDFSAVESLQ VLKLRWCKFE IPDSSPKGLK LLKTLSLMRT QVMVKTIDAI FNNCIHLESL ELIECRM DGI LSIRAQNHKK FKSLVVSFMP DLRHIRLDAP TLENYKYDGY VICVNILITN ALKEANLYYT RIRRLYHQKS DLVDTLRFYT RLTVLATTTI FLEALTKRYV GEGRENPPF KFENLTEFKI SFITPTFCTL FDIAEFLKEC PKLKQVVIDI QNFTFEPQMY FWEIHHKAQI QNTSNNNYLL KCLTDVKIIG YKGHWHELDI VEFFVKNAPS LKRLELQMPK NAKNDAHTPD VARIKLIKTI FSGVKVTEV
Source	Baculovirus
Target Names	At1g61320
Protein Names	Recommended name: FBD-associated F-box protein At1g61320
Expression Region	1-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.