



Recombinant Arabidopsis thaliana Probable germin-like protein subfamily 2 member 5 (At5g26700)

Product Code	CSB-EP524054DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O65252
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	EMLQDVCV ADLSNAVKVN GYTCKDSTQI TPEDFYFKGL ANIAATNTST GSVVTGANVE KLPGLNTLGL SMSRIDYAPN GLNPPHVHPR ASEIIFVLEG QLYVGFVTTA GKLIKLNK GDVFTFPKGL IHFQKNIANS PASVLAAFDS QLPGTQSLVA SLFGALPDDI LAKSFQLKHK QVKKIKLRYA PKK
Source	E.coli
Target Names	At5g26700
Protein Names	Recommended name: Probable germin-like protein subfamily 2 member 5
Expression Region	23-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.