



Recombinant Arabidopsis thaliana Putative pectate lyase 17 (At4g22090)

Product Code	CSB-MP524060DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O65457
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	AYNGYYGY SPAAPYPAE EPQNIMNPVD SCWRLKSDWD VNREDLADCA VGFGSSTLGG KKGNIYVVTN PYDNAQNPHP GSLRYGVIQA KPLWITFAKD MVITLANELM VNSYKTIDGR GAKVEIAYGP CITIQDVTNV IVHGSIHDC KPGKSGKQRS SPTHVGHKRG SDGDAITIFG SSNVWIDHCY LASCTDGLID VIHASTAITI SNNYFTQHDK VMLLGHNDNF VKDVKMKVTV AFNHFGPGLV ERMPVRRGY AHVANNRYDK WIMYAIGGSA DPTIFSEGN YFIASDKSYK EVTKREVKGG WNNWRWRTSN DVFKNGAFFV PSGYGSIPLP YSSAQRFTVA PGNLVPSLTA DAGPLNCNRN GPCY
Source	Mammalian cell
Target Names	At4g22090
Protein Names	Recommended name: Putative pectate lyase 17 EC= 4.2.2.2
Expression Region	23-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.