



Recombinant *Schizosaccharomyces pombe* Protein sof1 (sof1)

Product Code	CSB-EP524245SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74340
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MKVKITRGT SLTRLNDQDP VKRNLDPSLH PFERAREYTR ALNATKMDRM FAAPFLGQLG RGHQDGVYSL ARDTKTLIDC ASGSGDGAVK LWDASERCER WTSKAHEGIV RGLVFSNQQD VLSCASDRYV YMLNKQDGKV KRSYLGSSSL LDIDTSKGGD LFATSGENV S IWDYSRDTVP TKFEWGADTL PVVKFNYTET SVLASAGMDR SIVIYDLRTS SPLTKLITKL RTNSISWNPM EAFNFVAGSE DHNLYMYDMR NLKRALHVYK DHVSAVMSVD FSPTGQEFVS GSYDKTIRIY NVREGHSRDV YHTKRMQRVT AVKFSMDAQY IFSGSDDSNV RLWRARASSR ASIRSTREEN RLKYLDLRE RYKHIPEIRR IARHRHLPTN VKKAAEIKRE EINSLKRREE NIRRHSSKGA VPYEKERERH VVGIIQK
Source	E.coli
Target Names	sof1
Protein Names	Recommended name: Protein sof1 Alternative name(s): U3 small nucleolar RNA-associated protein sof1 Short name= U3 snoRNA-associated protein sof1
Expression Region	1-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.