



Recombinant *Schizosaccharomyces pombe* Uncharacterized protein C14G10.04 (SPCC14G10.04)

Product Code	CSB-EP524255SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74417
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSQVIQYSEA DLLYLSKSPL IKKPENLPEW ILPEAMKADR ERHIRQEKEM KRHDDDGRQY QSDRKFAKSK HDDIILGPPK LSFASSNENG RSVGRHDDL SLGNVYNGVA SLRYRNGASV SRSSSIGHSG STAPWSSVGR HNRKKDNEHR DEMEGLEKVM KQRAGNTGPL VANSTEEFEA WKRMMKASSG EAGAGNVITP GVTSTTGAPS GKASLSRAAS NSSTSARAGI SVDSLFGKRS AAQAMATVVP SVSTPGGTPP RFASPALVE SMVHAAPSTT TSSRFSKFFN GLAGPESKT STPNMSNPS PRVANTVVSP APIPTIGSAT IDKDAEGFQR VLAMLGRQAS STTGSPKVAL SQMPQVNVNP SNQDLASVKQ PSGFSPPSAV APPPGLGNHN FSNDSSFFR SLMTSDTRSV PPPPGFSTNA PKAVKSPEIS APPGLYRGLS SGASIPSAPP GFGYQQPSFP YSPGFPPSAY NQNRTGYGFP DTSRPPH
Source	E.coli
Target Names	SPCC14G10.04
Protein Names	Recommended name: Uncharacterized protein C14G10.04
Expression Region	1-497
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.