



Recombinant Schizosaccharomyces pombe 37S ribosomal protein mrp51, mitochondrial (mrp51)

Product Code	CSB-YP524256SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O42855
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	ASPKCLP ESNYSTSHLR FPRLQSITSK KNSRARGDWG LKRTLPKIKT RYICVRQQDS LYKQADFQSS AKFVRLVRNW FDAGFVRDPS EVALLQELLK ETGLKDTRDI NKVFSPSKQS DENDNFRRYG VSGGIQYSNV PLINSRVTPN CESNCSNIRR RLTAHISSD ASTAYYGLGG LILRLPRTVV FNRNRSSPLF SQEQRISYTV YKNGRLRLVP FSGPELETHL AYPSEDSEVE SFFNRDIQLV AKRDSSLDAP VSYQRGIMKF FMAPSQTSNN EANSSEKDNL SNS
Source	Yeast
Target Names	mrp51
Protein Names	Recommended name: 37S ribosomal protein mrp51, mitochondrial
Expression Region	14-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.