



Recombinant Schizosaccharomyces pombe Transcription factor nrm1 (nrm1)

Product Code	CSB-BP524262SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O42913
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDRSMEPLTP SRLNTLGERP TNEVYEGKG KNVQHLPIT PMQRPLGKEN AAPGTISPIA VRSRNVRAVE IADENACEEP VLKIKSVSST ESEEEKESST EIGEKQEKET HLEPKTPVQT NTNHHLDDI QCCAKNLRRLR LELAMYKVQV NQTFSPQLDL PIVAKTKLHN CPNSEPVTSI WNQRSLSSGK PPSLHLSGNR RLSMGSPTKS IYDQNGLTTP RPIGSDDLTH MYDPYTSPLR TPSRTLRSR SHYLWVRHGK LTRSVSLLQH KTPRRIRPKS LSKSNSTPLK HLSAQKPNNS YYTGPPTPVS ISNTPENIHP SSSEVRRIAS HSKQFSDYGL IR
Source	Baculovirus
Target Names	nrm1
Protein Names	Recommended name: Transcription factor nrm1 Alternative name(s): Negative regulator of MBF targets 1
Expression Region	1-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.