



Recombinant Pyruvate dehydrogenase E1 component subunit beta, mitochondrial (C04C3.3)

Product Code	CSB-YP524287CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O44451
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	ASTMTVRDA LNQAMDEEIK RDDRVLMEGE EVAQYDGAYK ISKGLWKKHG DKRVVDTPIT EMGFAGIAVG AAFAGLRPIC EFMTFNFSMQ AIDQIINSAA KTYYSAGRV PVPIVFRGPN GAAAGVAAQH SQDFSAWYAH CPGLKVVCPY SAEDAKGLLK AAIRDDNPVV FLENEILYGG SFPVGDEVLS DDFVVPYIGKA KIERAGDHVT IVSYSRGVEF SLEAAKQLEA IGVSAEVINL RSLRPFDFES IRQSVHKTHH LSVVETGWPF AGIGSEIAAQ VMESDVFQDL DAPLLRVTGV DVPMPYQTL EAAALPTAEH VVKAVKKS LN IA
Source	Yeast
Target Names	pdhb-1
Protein Names	Recommended name: Pyruvate dehydrogenase E1 component subunit beta, mitochondrial Short name= PDHE1-B EC= 1.2.4.1
Expression Region	22-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.