



Recombinant Ephrin-4 (efn-4)

Product Code	CSB-MP524288CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O44516
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	DEHIVYWNST NSLFRNRQPT IEVRMGDVVR FVCPDNEEGR NDGEYLIVYE VTEFAMDDCA LESHGREVIR CAPEGTAEKV LRTQQLSGGR REDWKKQKVP PKNVAQLIRQ LNPIPNGKEY QPGQTYYYMT TSTGKANGTN HRMYGLCESQ NMRLSMKVSA SQPHPTRRAP TRRQEDFVTT ASAEIMGGQE DESDNDNAH LLPRDLEGST NPKFRRPSQL ETAGVENQQF MKVVQMAQAG KTGTFENEKE AIAQKSSEKD GWHPVNVQYV ADLMNNAYQN ADERISYQRD FEIHEENDLA VKSLEYSSS
Source	Mammalian cell
Target Names	efn-4
Protein Names	Recommended name: Ephrin-4 Alternative name(s): Protein male abnormal 26
Expression Region	21-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.