



# Recombinant Arabidopsis thaliana Protein disulfide-isomerase 2-3 (PDIL2-3)

<b>Product Code</b>	CSB-MP524331DOA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O48773
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	LYGSSS PVVQLTASNFKSKVLNSNGV VLVEFFAPWC GHCKALTPTW EKVANILKGV ATVAIDADA HQSAAQDYGI KGFPTIKVFV PGKAPIDYQG ARDAKSANF AYKQIKGLLS DRLEGKSKPT GGSKEKKSE PSASVELNAS NFDDLVIEN ELWIVEFFAP WCGHCKKLAP EWKRAAKNLQ GKVKLGHVNC DVEQSIMSRF KVQGFPTILV FGPDKSSPYP YEGARSASAI ESFASVELS SAGPVEVTEL TGPDMMEKKC GSAAICFISF LPDILDSKAE GRNKYLEMLL SVAEKFKKQP YSFMWVAVT QMDLEKRVNV GGYGYPAMVA MNVKKG VYAP LKSAFELQHL LEFVKDAGTG GKG NVP MNGT PEIVKTK EWD GKDGELIEED EFSLDELMGG DDAVGSKDEL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PDIL2-3
<b>Protein Names</b>	Recommended name: Protein disulfide-isomerase 2-3 Short name= AtPDIL2-3 EC= 5.3.4.1 Alternative name(s): Protein disulfide-isomerase 5-2 Short name= AtPDIL5-2 Protein disulfide-isomerase 9 Short name= PDI9
<b>Expression Region</b>	25-440
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.