



# Recombinant Arabidopsis thaliana Gibberellin 2-beta-dioxygenase 8 (GA2OX7)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP524350DOA   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | O49561  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Arabidopsis thaliana (Mouse-ear cress)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MDPPFNEIYN NLLYNQITKK DNDVSEIPFS FSVTAVVEEV ELPVIDVSRL<br>IDGAEEREK CKEAIARASR EWGFFQVINH GISMDVLEKM RQEQRVFRE<br>PFDKKSSEK FSAGSYRWGT PSATSIRQLS WSEAFHVPMT DISDNKDFTT<br>LSSTMEKFAS ESEALAYMLA EVLAEKSGQN SSFFKENCVR NTCYLRMNR<br>PPCPKPSEVY GLMPHTDSDF LTILYQDQVG GLQLIKDNRW IAVKPNPKAL<br>IINIGDLFQA WSNMGYKSVE HRVMTNPKVE RFSTAYFMCP SYDAVIECSS<br>DRPAYRNFSF REFRQQVQED VKKFGFKVGL PRFLNHVY           |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | GA2OX8  |
| <b>Protein Names</b>     | Recommended name: Gibberellin 2-beta-dioxygenase 8 EC= 1.14.11.13<br>Alternative name(s): GA 2-oxidase 8 Gibberellin 2-beta-hydroxylase 8<br>Gibberellin 2-oxidase 8  |
| <b>Expression Region</b> | 1-338   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |