



Recombinant *Zygosaccharomyces rouxii* Mitogen-activated protein kinase HOG1 (HOG1)

Product Code	CSB-YP524579ZBK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O59854
Product Type	Recombinant Protein
Immunogen Species	<i>Zygosaccharomyces rouxii</i> (<i>Candida mogii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MATHEEFIRT QVFGTVFEIT NRYTDLNPVG MGAFLVCSA TDTLAGQPVA IKKIMKPFST AVLAKRTYRE LKLLKHLRHE NLICLQDIFL SPLEDIYFVT ELQGTDLHRL LQTRPLEKQF VQYFLYQILR GLKYVHSAGV IHRDLKPSNI LINENCDLKI CDFGLARIQD PQMTGYVSTR YYRAPEIMLT WQKYDVEVDI WSAGCIFSEM IEGKPLFPGK VHVHQFSIIT DLLGSPPRDV IITICSEDTL KFVTSPLHRD PVPFQERFKA VEPDAVDLLG RMLVFDPKKR ITAADALVHP YLAPYHDPTD EPIAEAQFDW DFNDADLPVD TWRVMMYSEI LDFHKIGGSD GQIDTNAAFD DQVAAATAGP MLPLWLNSNI NINNSHPVNI QVQLLRHQQP RITVVYQ
Source	Yeast
Target Names	HOG1
Protein Names	Recommended name: Mitogen-activated protein kinase HOG1 Short name= MAP kinase HOG1 EC= 2.7.11.24 Alternative name(s): ZrHOG1
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.