



Recombinant Schizosaccharomyces pombe Uncharacterized mannosyltransferase C19C7.12c (SPBC19C7.12c)

Product Code	CSB-EP524595SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O60160
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVRLPRKFKR VLLLVLTL VVFVRFKKQY IPTISVFEGS LIDNRDTLSY FNISNLEPSE RSEWLPNKRV NAAFVTLARN EDLNDLLKSI RKLEKTFNHK YHYGWVFLNN EEFSDDEFKEH VIEAVSGKCE FGLVPQEQWS IPEYIDQDKM HENWKKLQEL GILYADKESY RHMCRFESGF FYRHPLVQKY EYYWRVEPSV DFFCDLDFDP FAYMKENKA YAFTITVTEY SETIPSLWPS TKEFIKMHPN ALHPNNALNF ISNDDGETYN GCHFWTNFEI AKVDFWESEV YSKYFDYLDK SGNFFYERWG DAPVHSIAVS LFADRDNIHF FNEIGYWHPG SGHCPLDDAT RAKDCDPYE SIDYNGWSCL DKFYTAFEMP FPENWSFYSH
Source	E.coli
Target Names	omh1
Protein Names	Recommended name: Uncharacterized mannosyltransferase C19C7.12c EC=2.4.1.-
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.