



Recombinant Arabidopsis thaliana Probable inactive poly [ADP-ribose] polymerase SRO3 (SRO3)

Product Code	CSB-BP524629DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O64592
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAAQVEIEDQ ESVTNLDNGE IINPISDNAP NFSGDATILL REATFEHNLI KNCFLSGMGS FATETTIVTV RKILTQRLIT TKAKFAVFKL FTEAMKRKNN GYANIRYGWY SGSKEEIDRV ITYGFSNREI KKVENDVGSH GVGIHLVHHR YSLAAALVGE GDEEGIKNIL LCRVILGKPE QIVTGSKQSY PSSNQFDSGV DNLENPRKYV IWSCNMNSYI LPTYIVSFKS HLLRGLIGRA RSPCVSFSVL MSILSKSLDA ARMNILITSY DDFRKRKLRR EQLVRKIREV VGDNLLFKIL KNQRR
Source	Baculovirus
Target Names	SRO3
Protein Names	Recommended name: Probable inactive poly [ADP-ribose] polymerase SRO3 Alternative name(s): Protein SIMILAR TO RCD ONE 3
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.