



Recombinant Arabidopsis thaliana Coatomer subunit epsilon-2 (At2g34840)

Product Code	CSB-EP524633DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O64748
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSAMGAGPDH LFNLRNMFYL GAYQTAINNS EIANLSPENA VERDCLVFRS YIALGSYQLV ISEIDESAAAT PLQAVKLLAM YLSTPQNKES TISSLKEWLA DSTIGNNDTL RLIAGIIFMH EEDYNETLKH THAGGTMDLY ALNVQIFIKM HRAEYAEKQL RVMQQIDEDH TLTQLASAWL NLAVGGSKIQ EAYLIFEDFS EKYPMTCLIL NGKAVCCMQM GNFDEAETLL LEALNKDAKD PETLANLVVC SLHVGKSSSR HLSQLKLSHP EHILVKRVSS AEDNFERAVQ LVA
Source	E.coli
Target Names	At2g34840
Protein Names	Recommended name: Coatomer subunit epsilon-2 Alternative name(s): Epsilon-coat protein 2 Short name= Epsilon-COP 2
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.