



Recombinant Arabidopsis thaliana Putative ubiquitin-like-specific protease 1B (ULP1B)

Product Code	CSB-EP524649DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O65278
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MALVPKIRNR TISDKPVSSS CTRTSSFRQR LLASCKRFLR FKLLVRASSV SLKKKAARCS QSQQVVKHND DRIEEHSSSN RGFKDGNNFL VSSTETKKLV AFCGEPFQPL NEEVALVNS ALSKRNRKKI LVSHKNSNID ISGETLQCLR PNQWLNDVDT NLYLELLKER QTRDPQKYFK CHFFNTFFYV KLVSGSGYNY KAVSRWTTKR KLGYDLIDCD IIFVPIHIDI HWTLGVINNR ERKFVYLDL FTGVGHTILN AMAKYLVEV KQKSQKNIDV SSWGMEYVEE RPQQQNGYDC GMFMLKYIDF YSRGLSLQFS QKDMPYFRLR TAKEILRLRA D
Source	E.coli
Target Names	ULP1B
Protein Names	Recommended name: Putative ubiquitin-like-specific protease 1B EC= 3.4.22.68
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.