



Recombinant Chicken Glutamine synthetase (GLUL)

Product Code	CSB-BP009553CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P16580
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	MATSASSHLS KAIKHYMYMKL PQGEKVQAMY IWIDGTGEHL RCKTRTL DHE PKSLEDLPEW NFDGSSTFQA EGSNSDMYLR PAAMFRDPFR KDPNKLVLCE VFKYNRQSAD TNLRHTCRRI MDMVSNQHPW FGMEQEYTLT GTDGHPFGWP SNCFPGPQGP YYCGVGADKA YGRDIVEAHY RACLYAGVKI GGTNAEVMPA QWEFQVGPCE GIEMGDHLWI ARFILHRVCE DFGVIVSFPD KPIPGNWNGA GCHTNFSTKN MREDGGLKHI EEAIEKLSKR HQYHIRAYDP KGGLDNARRL TGFHETSSIH EFSAGVANRG ASIRIPRNVG HEKKG YFEDR GPSANCDPYA VTEALVRTCL LNETGDEPFE YKN
Source	Baculovirus
Target Names	GLUL
Protein Names	Recommended name: Glutamine synthetase Short name= GS EC= 6.3.1.2 Alternative name(s): Glutamate decarboxylase EC= 4.1.1.15 Glutamate-- ammonia ligase p42
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.