



Recombinant Schizosaccharomyces pombe Uncharacterized protein C162.12 (SPCC162.12, SPCC1753.06c)

Product Code	CSB-EP524854SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74428
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MERPSLSRRT SSSTVSTDGE GVYSRSTKER KRNFI VNSRL KRGGGHVRVN RSGRLGSVTM RPSALTRAHS QNPNSSLVNN SSAVSHLTKQ RSLNDLHELN GPKHVAKNGL IPLTQRKPNC VWDDAPVDND STAGNLDSDS ALPTPSVTTN EAADSSRASS PVTRVAVHD NKKKIINSNI SNAPPFNNTD VQASARPPAA GQDDSAADAS TTKSSPVHNE VMAEPLPHSN NREVTQATNQ PKWQIHSGSD IASEPPTLSR GNSL SLLANR KVPSTNVKKS QAELYDLGSS TSRTQQKLLI QRASSKFDIV EDDMDTNPSK RFSNPHTKHI MDLVRTQYRN VLRTRELIPE FLEKIRSSYS NNQNFDSQNA FNTSAAGTAG TREETISNGQ NGIVASAETS KKDDGVQSAS LNASMSARSH ARQRSIHVPK TRKDTDYESI HQKLLQLWSQ G
Source	E.coli
Target Names	tco89
Protein Names	Recommended name: Uncharacterized protein C162.12
Expression Region	1-451
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.