



# Recombinant *Schizosaccharomyces pombe* tRNA (guanine-N (7)-)-methyltransferase subunit trm82

<b>Product Code</b>	CSB-YP524900SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O74863
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSEQRVFKHP CQFLTWNSKH NYIVCCSGPY LLGFSCSTGE KIFEHCYRDN INEKHREAAA YGEAIRQVAF SKDYSRMATV SEDKCLRLWD STQPKIELL YQKNIPKRCA DLCFAGSNEI VFGDKFGDVY CVDENWFTTS EVTEEKKSNV VEGKQEPVNN DTLKDSKLQK LEPIMGHVSI LTQLIVAQNP QNSKEEIIIT SDKDEHIRIS RFPNAFVIEG FCLGHEDFVS RMSLYDNRTL ISGGGDNHVF VWDLENFKCL DAFDLRSAFS TYLSLNQPMV VSVILPIFKR QLVAFACEGM AGLIFAKVTP EKRLLFHSAL KLSGPVLDVAV LLDTDTDQIL ISLDSSFTFG ACFECVKFDE GNAAVLTKPD VIKRIESDGL ISTEKPFPCPL AQIHTLRKNH SKFIRSVETG TSSPSVESKD N
<b>Source</b>	Yeast
<b>Target Names</b>	trm82
<b>Protein Names</b>	Recommended name: tRNA (guanine-N(7)-)-methyltransferase subunit trm82 Alternative name(s): Transfer RNA methyltransferase 82
<b>Expression Region</b>	1-421
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.