



Recombinant Schizosaccharomyces pombe 26S protease regulatory subunit 6B homolog (rpt3)

Product Code	CSB-EP524905SXV-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O74894
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MNSTEEIDLY NRWKALERQL EMLDLQEGFI KEDCKSLKRE LIRAQEEVKR IQSVPLVIGQ FLEAIDQNTA IVGSTTGSNY VVRILSTLDR ELLKPSASVA LQRHSNALVD ILPPEADSSI SMLRPDERPD VSYADVGLD VQKQEVREAV ELPLTQGDLY RQIGIDPPRG VLLYGPPGTG KTMLVKAVAN STAA NFIRVV GSEFVQKYL G EGPRMVRDVF RMARENAPAI IFIDEIDAIA TKRFDAQTGA DREVQRILIE LLTQMDGFDQ GANVKVIMAT NRADTLDPAL LRPGR LDRKI EFPSYRDRRQ RRLVFQTITA KMLLSPEVDL DTFIMRPDAS SGAQIAAIMQ DAGLLAVRKS RGVILQSDIE EAYARAVKPT DMEQFAFYK
Source	E.coli
Target Names	rpt3
Protein Names	Recommended name: 26S protease regulatory subunit 6B homolog
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.