



Recombinant Arabidopsis thaliana GDSL esterase/lipase At2g38180 (At2g38180)

Product Code	CSB-BP524966DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O80443
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	DRGWGATL ADLYSRTADI ILRGYAGWNS RFALKVLHQV FPKDAVIQPS LVIVYFGGND STHPHPSGHG PHVPLSEFIE NMRKIGEHL SLSDKTRVIF LTPPPMNEKQ IEIVFGDAIK GRSNELCRPY AEELNLCRE INVKGIDIWT AIQQQDDWLN SCFTDGIHFT AKASEIVVKE ILKVLRGADW KPSLYWKSLP VEFPFDFDAP NSISLHDLEL TRNNHFESPH LVSLCEQELT RNEQLEPPHP VSLCDHELTR NEQLEPPHPV SLCDHELTQN EQLEPPQPTA RL
Source	Baculovirus
Target Names	At2g38180
Protein Names	Recommended name: GDSL esterase/lipase At2g38180 EC= 3.1.1.- Alternative name(s): Extracellular lipase At2g38180
Expression Region	23-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.