



# Recombinant Acetobacter pasteurianus Type-2 restriction enzyme ApaLI (apaLIR)

<b>Product Code</b>	CSB-MP525151AVW
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O52703
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Acetobacter pasteurianus (Acetobacter turbidans)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTTRQRLSAE RSQQLTRLLT ITKTANMRAL MEASELAKVI ALVAVDIGKS DEMARAFPVL WPKISPQQEY YATAVDWFTN PDETVTSTFDV VDMLDAGTSL DQDFMTYLKC LTELHKRRRK YGLILQRQPL PTMVQVSPRA LMEYGPDFPP EALASWLTWR KFFYDLNRS AQETGYLFEP ILAAAIGGEA KSARERVVRR TDDPTKGRQV DCWKVLPDGT PLAYELKLRV TIAASGQGRF GEELSFARDC SSSGAKPILV VLDPTENDKL TGLQAAYREV GGAAYVGDAW WAHLEDEAGA TMAFSIERYV RVPVASVSSF ERVIEGDATK RSLILQDLQA RLDGNELTIS LGGHQRLVER HEDQSLAADG DDDSE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	apaLIR
<b>Protein Names</b>	Recommended name: Type-2 restriction enzyme ApaLI Short name= R.ApaLI EC= 3.1.21.4 Alternative name(s): Endonuclease ApaLI Type II restriction enzyme ApaLI
<b>Expression Region</b>	1-375
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.