



Recombinant Schizosaccharomyces pombe Chitin synthase regulatory factor 1 (chr1)

Product Code	CSB-BP525282SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O59732
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MPASTSIPET PSKLPDALLV NSDPASKVDN VSVKIEQESP LALEGSPLPT VSDVLPETNN AHEEVSSSSW SLSSSDSIMS YDSLSDDVSV LSFLDCPLTT FETAPSIASF SNSAADVSSH SRSSSSWSAK LSRFTKHASK PSLSSNSSDS SFSKSGSESN ALSSRGSFEL EPHGIHRHVS RAKRLLSKFY SKFHHKKEDS SFDPLAPLVF APNTSRVLRV TNEANAS AIS LEKATSLSSQ EQPGSEIKKE LSLYDHEFDQ ILDLAIKAKK DGNLENAIEF LEYIMPEVAS QQNFVPYELA ELYKQRGTSQ DLKSILPLYM LAASLGHDRS SFLVGEAFFY GTYGARENKL RALQYYHLAN DKGNAAML LCKLYLRGLP GHIFPSSRRA FEYAHRAAML GHAPACYVLG KFYETGVGCV KDLAKSEAGF RAGLINDSPI TDADALQIAS TLVLLQ
Source	Baculovirus
Target Names	chr1
Protein Names	Recommended name: Chitin synthase regulatory factor 1 Alternative name(s): Chs four homolog 4
Expression Region	1-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.