



Recombinant Schizosaccharomyces pombe Zinc finger protein C550.15c (SPCC550.15c)

Product Code	CSB-BP525290SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O59811
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSTSFAC TTC TVAFNNAESQ KIH WKSDWHH YNLKRKVASL PPLSAEVFAG KILSIQKQNE EVQKKAEFYQ NCEVCNKKFY SEGAYSSHMA SKKHRDNL SK FQRNSRIKKL QSEDASSIAS STLSMGEPVV DSEIEEEEEDL ASQLTSRAIS LSNLSLHGREG SEPSKTELAT SIPQSNEASK SHLFTQEPTP EEIEAELARR SSQRLSPRDC LFCAASFSSF DTCKKHKMKAS HSLYIPEREY LVDEPSLFDY LAEKISIGFT CLTCNREFKS LEAVRAHMQQ KGHTSIAYDT EDEQLELSDF YDFTTSYPDY AVKQDETVE EDGSSGEGDW EDVSDDSDNS SLDSLEMGRV PIADEYELHL PSGNRVGHRS LSR YFRQNLH SSSTAVGDGA SIHQNVARRA MSGNARAYRQ AVETSIAGVR DGRKNYSASH IKS FQDQRRR EEFANKMGIK NNTKKHFRDA LLQ
Source	Baculovirus
Target Names	SPCC550.15c
Protein Names	Recommended name: Zinc finger protein C550.15c
Expression Region	1-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.