



# Recombinant Schizosaccharomyces pombe Carbamoyl-phosphate synthase arginine-specific small chain (arg5)

<b>Product Code</b>	CSB-MP525302SXV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O60060
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLRFLKPFPL RFGKRFYSKV PPVTNHERIL PKQPSFPTAP AQNEIATLTI RNGPIFHGTS FGANRNVSGE AVFTTSPVGY VESLTDPSYK QQILIFTQPL IGNYGVDPCK KRDENGLLRH FESPHIQCAG VVVNDYATKY SHWTAVESLG EWCAREGVAA ITGVDTRAIV TFLREQGSSL AKISIGEEYD ANDDEAFINP EEVNLVSQVS TREPFFVSGG DGMLNIAVID CGVKENILRS LVSARGASVTV PPFDYPIQNV ASNYDGIFLT NGPGDPHTLT KTVNNLRELM NTYNGPIMGI CMGHQLLALS TGAKTIKLY GNRGHNIPAL DIASGNCHIT SQNHGYAVDA STLPAEWKAT WTNLNDQSNE GIAHVS RPIS SVQFHPEAR GPMDFYLF DYI KEAIKYQ KSR TA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	arg5
<b>Protein Names</b>	Recommended name: Carbamoyl-phosphate synthase arginine-specific small chain Short name= CPS-A EC= 6.3.5.5 Alternative name(s): Arginine-specific carbamoyl-phosphate synthetase, glutamine chain
<b>Expression Region</b>	1-415
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.