



Recombinant Schizosaccharomyces pombe 54S ribosomal protein L23, mitochondrial (mrpl23)

Product Code	CSB-MP525306SXV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O60118
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	A SGIATTLMGK HKPIYHPAAD CGDVVVVTDG SEIGISGRKL ENHKYYSHSG QPGHLKEWTM EEMAAKRGHK DLLRRAVGGM LPRNKLWDRR MKRLYIYDGA EHPYKANIFR SYHNPVSSQI AKELFSK
Source	Mammalian cell
Target Names	mrpl23
Protein Names	Recommended name: 54S ribosomal protein L23, mitochondrial
Expression Region	30-157
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.