



Recombinant Probable selenide, water dikinase (seld-1)

Product Code	CSB-EP525328CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O62461
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MNRIERILEG FDPVSNGLDE DFVLTKLTGM KGCGCKVPRN VLLQLLQTFK TDLVINNDEV DIGLDSCVIP LRHPGLRLVQ TTDDFFYPLID DPYIMGRVTC ANVLSDLYAM GVSECDNMLM LLAVAIDLNE KQRDIVVPLF IQGFKDAADE AGTKIRGGQT VRCPWLLLGG VATSVAHESE IIKVDQAVPG DVLILTKPIG GQVAVNSYEW IKKKNKIEE LNLEIPKIEK AFKQVCEQMS RLNRNAAKLL HKYDAHSSTD VTGFLLGHA ENLARVQKQP MEFIIEKLPI IEYMDEIADK MIAKGGEGFK LYQGTAETS GGLLIAMSEE NAKKYIAELS SLDNAPAWII GKVTAKTTDS SIARILPDAV RISVPSHI
Source	E.coli
Target Names	seld-1
Protein Names	Recommended name: Probable selenide, water dikinase EC= 2.7.9.3 Alternative name(s): Selenium donor protein Selenophosphate synthase 1
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.