



Recombinant Arabidopsis thaliana Protein TIFY 5B (TIFY 5B)

Product Code	CSB-BP525358DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O64687
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MIIIIKNCDK PLLNFKEMEM QTKCDLELRL LTSSYDSDFH SSLDESSSSE ISQPKQESQI LTIFYNGHMC VSSDLTHLEA NAILSLASRD VEEKSLSLRS SDGSDPPTIP NNSTRFHYQK ASMKRSLHSF LQKRSLRIQA TSPYHRYR
Source	Baculovirus
Target Names	TIFY
Protein Names	Recommended name: Protein TIFY 5B Alternative name(s): Jasmonate ZIM domain-containing protein 7
Expression Region	1-148
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.