



Recombinant Arabidopsis thaliana Putative FBD-associated F-box protein At1g61330 (At1g61330)

Product Code	CSB-BP525363DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O64787
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MSEASNKQMK LIKRLSDELV ECILSFLPVQ STLQHRVLSK RYRDTWKLSR DLDFSGIFSV RRSQSEAVRI IEDVFIHQHG PEIDRFVLSL DHTGVEDKIL SWVDTCLRKN IKELVLDLDFSS SKKVMEIPLN FSSSETLTVL KLQWCRFEIP NNLPKGLRLL KTLSLMRTEV TNEMIDSIFN NCIHLESLEL VKCQMSGNGI LTIYAHDHKK FKSLVVSCMP NLLSIFLDAP TLECYKYDGY AKTINVLKVN ALEEAEFHYS RSKRRHSSDT VVDSMLACSA VHALTTTNil LEAITRYLK GKLEKPLFKF ENLREFKIFF KAPTFCTLFD IAEFLKECPK LEQFIIDIQN FTFEPQMYFW EIIHKAQIQN TSNNNYLLKC LTDVKIIGYK GHWHELDIVE FFVKNAPSLK RLELQMPKNA KND AHTPDVA RIKLIKTIFS GVKVTEV
Source	Baculovirus
Target Names	At1g61330
Protein Names	Recommended name: Putative FBD-associated F-box protein At1g61330
Expression Region	1-447
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.