



# Recombinant Arabidopsis thaliana Putative FBD-associated F-box protein At1g61330 (At1g61330)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP525363DOA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                   |
| <b>Uniprot No.</b>       | O64787                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Immunogen Species</b> | Arabidopsis thaliana (Mouse-ear cress)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Sequence</b>          | MSEASNQMK LIKRLSDELV ECILSFLPVQ STLQHRVLSK RYRDTWKLSR<br>DLDFSGIFSV RRSQSEAVRI IEDVFIQHQG PEIDRFVLSL DHTGVEDKIL<br>SWVDTCLRKN IKELVLD FSS SKKVMEIPLN FSSSETLTVL KLQWCRFEIP<br>NNLPKGLRLL KTLSLMRTEV TNEMIDSIFN NCIHLELEL VKCQMSGNGI<br>LTIYAHDHKK FKSLVVSCMP NLLSIFLDAP TLECYKYDGY AKTINV LKVN<br>ALEEAEFHYS RSKRRHSSDT VVDSMLACSA VHALTTT NIL LEAITRYLK<br>GKLEKPLFKF ENLREFKIFF KAPTFCTLFD IAEFLKECPK LEQFIIDIQN<br>FTFEPQMYFW EIIHKAQIQN TSNNNYLLKC LTDVKIIGYK GHWHELDIVE<br>FFVKNAPSLK RLELQMPKNA KND AHTPDVA RIKLIKTIFS GVKVTEV |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Target Names</b>      | At1g61330                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Protein Names</b>     | Recommended name: Putative FBD-associated F-box protein At1g61330                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Expression Region</b> | 1-447                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                          |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                   |